



Laser Diode Driver

LDD605



Revision 2.11

Limitation of Liability

MOG Laboratories Pty Ltd (MOGLabs) does not assume any liability arising out of the use of the information contained within this manual. This document may contain or reference information and products protected by copyrights or patents and does not convey any license under the patent rights of MOGLabs, nor the rights of others. MOGLabs will not be liable for any defect in hardware or software or loss or inadequacy of data of any kind, or for any direct, indirect, incidental, or consequential damages in connections with or arising out of the performance or use of any of its products. The foregoing limitation of liability shall be equally applicable to any service provided by MOGLabs.

Copyright

Copyright © MOG Laboratories Pty Ltd (MOGLabs) 2016 – 2025. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form or by any means, electronic, mechanical, photocopying or otherwise, without the prior written permission of MOGLabs.

Contact

For further information, please contact:

MOG Laboratories P/L
49 University St
Carlton VIC 3053
AUSTRALIA
+61 3 9939 0677
info@moglabs.com
www.moglabs.com

Santec LIS Corporation
5823 Ohkusa-Nenjozaka, Komaki
Aichi 485-0802
JAPAN
+81 568 79 3535
www.santec.com

Preface

The MOGLabs LDD Laser Diode Driver is a compact high-power laser diode driver. It provides up to 8 A diode injection current with low noise, a temperature controller, a 60 W Peltier TEC driver, and a 160 V piezo driver. It can be operated via front-panel controls, or using a computer communications interface (TCP/IP or USB) with simple text-based commands. It has been designed for tapered amplifier and broad area laser diodes commonly used in atomic physics research laboratories, where stability and low noise are essential.

We hope that the LDD meets and exceeds your expectations. Please let us know if you have any suggestions for improvement in the LDD or in this document, and check our website from time to time for updated information.

MOGLabs www.moglabs.com

Safety Precautions

Safe and effective use of this product is very important. Please read the following safety information before attempting to operate your laser. Also please note several specific and unusual cautionary notes before using the MOGLabs LDD, in addition to the safety precautions that are standard for any electronic equipment or for laser-related instrumentation.

CAUTION – USE OF CONTROLS OR ADJUSTMENTS OR PERFORMANCE OF PROCEDURES OTHER THAN THOSE SPECIFIED HEREIN MAY RESULT IN HAZARDOUS RADIATION EXPOSURE

Laser output can be dangerous. Please ensure that you implement the appropriate hazard minimisations for your environment, such as laser safety goggles, beam blocks, and door interlocks. MOGLabs takes no responsibility for safe configuration and use of your laser. Please:

- Avoid direct exposure to the beam.
- Avoid looking directly into the beam.
- Note the safety labels and heed their warnings.
- When the laser is switched on, there will be a short delay of two seconds before the emission of laser radiation, mandated by European laser safety regulations (IEC 60825-1).
- The STANDBY/RUN keyswitch must be turned to RUN before the laser can be switched on. The laser will not operate if the keyswitch is in the STANDBY position. The key cannot be

removed from the controller when it is in the clockwise (RUN) position.

- To completely shut off power to the unit, turn the keyswitch anti-clockwise (STANDBY position), switch the mains power switch at rear of unit to OFF, and unplug the unit.
- When the STANDBY/RUN keyswitch is on STANDBY, there cannot be power to the laser diode, but power can still be supplied to the laser head for temperature control.

CAUTION The supply must include a good ground connection. The laser chassis connected to the controller should be in good electrical contact to an optical table or other surface, which in turn should be connected to the mains power supply electrical ground.

CAUTION To ensure correct cooling airflow, the unit should not be operated with cover removed.

WARNING The internal circuit boards and many of the mounted components are at high voltage, with exposed conductors, in particular mains supply to various sections of the power supply. The unit should not be operated with cover removed.

NOTE The MOGLabs LDD is designed for use in scientific research laboratories. It should not be used for consumer or medical applications.

Protection Features

The MOGLabs LDD includes a number of features to protect you and your laser.

- Softstart** A time delay (3 s) followed by linearly ramping the diode current (1 s/A max).
- Circuit shutdown** Many areas of the circuitry are powered down when not in use. The diode current supplies may be without power when the unit is in standby mode, if an interlock is open, or a fault condition is detected.
- Current limit** Sets a maximum possible diode injection current.
- Cable continuity** If the laser is disconnected, the system will switch to standby and disable all laser supplies. If the laser diode, TEC or temperature sensor fail and become open-circuit, they will be disabled accordingly.
- Short circuit** If the laser diode, TEC or temperature sensor fail and become short-circuit, they will be disabled accordingly.
- Temperature** If the detected temperature is outside normal operating temperature, the temperature controller is disabled.
- Internal supplies** If one of the internal DC power supplies (+3.3, +5, ± 14 V) is 1 V or more below its nominal value, the respective components are disabled.
- Protection relay** When the power is off, or if the laser is off, the laser diode is shorted via a normally-closed solid-state relay at the laser head.

- Emission indicator** The MOGLabs controller will illuminate the current warning indicator LED immediately when the laser is switched on. There will then be a delay of at least 3 seconds before actual laser emission.
- Mains filter** Protection against mains voltage supply transient peaks.
- Key-operated** The laser cannot be powered unless the key-operated STANDBY switch is in the RUN position, to enable protection against unauthorised or accidental use. The key cannot be removed from the controller when it is in the clockwise (RUN) position.
- Interlocks** The main unit has an external interlock to allow the laser to be disabled via a remote switch, and a laser head cover interlock can be connected to disable laser output if the cover is removed.
- Seed input** A photodetector in the laser head can be used to detect the presence of a seed laser.

Contents

Preface	i
Safety Precautions	iii
Protection Features	v
1 Introduction	1
1.1 Front panel controls	1
1.2 Rear panel connections	2
1.3 Internal DIP switches	4
2 Menus	5
2.1 Main menu	5
2.2 Settings menu	6
2.3 TEC menu	6
2.4 Laser menu	7
2.5 Ramp settings	8
2.6 Interlock	9
2.7 Ethernet settings	9
2.8 Display	10
2.9 System information	10
3 Operation	13
3.1 Getting started	13
3.2 Interlocks	14
4 MOGLDD application	17
4.1 Device discovery	17
4.2 Main display	18
4.3 Menu options	19

4.4	Firmware update	22
A	Specifications	25
B	Troubleshooting	29
B.1	TEMP indicator	29
B.2	CURRENT indicator	29
B.3	Error states	30
C	Laser head board	33
C.1	B1055/B1056 headboard	33
D	Connector pinouts	37
D.1	Current	37
D.2	Temperature	38
D.3	Interlock	39
D.4	Piezo	39
E	Command language	41
E.1	General functions	41
E.2	Temperature control settings	43
E.3	Current control settings	45
E.4	Headboard settings	45
E.5	Piezo driver settings	46
E.6	Ethernet settings	47
F	Communications	49
F.1	Protocol	49
F.2	TCP/IP	50
F.3	USB	50


1. Introduction

The MOGLabs LDD is a compact high-power laser diode driver, providing low-noise injection current for the diode, a temperature controller with Peltier TEC output, and a high voltage piezo driver.

Connections between the LDD and laser are split across three cables: one for temperature control and power, one for diode current, and a third for the piezo. The connections are designed to interface directly with MOGLabs MOA and MSA optical amplifiers. See Appendix C for interfacing the LDD with non-MOGLabs laser products.

The device can be controlled via the front-panel controls and display (see below) or remotely via simple text commands communicated through USB or ethernet. Please refer to appendix F for further information on setting up computer communications, and appendix E for details on the control commands and their syntax.

1.1 Front panel controls



STANDBY/RUN

In STANDBY mode, both the temperature controller and diode current are disabled.

In RUN mode, the temperature controller is placed into standby, ready to be activated. On first power-up and during the boot pro-

cedure the LDD will halt if the key is already in the RUN position, to ensure that the laser cannot be accidentally powered if there is a power failure. The keyswitch should be set to STANDBY and then to RUN to initiate user control. Place the unit back into STANDBY if further operation is not desired.

TEMP OFF/ON Temperature controller enable. The STANDBY/RUN key switch must first be in the RUN position. If the unit fails to start the temperature controller on (indicator does not show green), see Appendix B.

CURRENT OFF/ON Diode current enable. The STANDBY/RUN key switch must first be in the RUN position. The temperature controller must also be on and operating.

If the current controller fails to switch to STANDBY mode (indicator does not show yellow), see Appendix B.

ADJUST The ADJUST rotary encoder allows variation of any editable value on the display. The knob can be pressed to step through different digits of the value being edited, to allow for fine and coarse control.


←, → The two buttons below the adjust knob allow selection of the digit of a value being edited.

←, ↑, ↓, OK Menu control buttons. These are used to navigate the menu system allowing users to go back one menu level (←), up and down between lines on any given menu (↑ or ↓), and to move forward in the menu system or to run system functions (OK).

1.2 Rear panel connections

IEC power in The unit supports a wide range of input voltages and frequencies from 80-264VAC 47-63Hz.

Fan The fan speed is temperature-controlled.



RAMP OUT	SMA connector. Represents the ramp applied to the piezo, when that functionality is enabled.
CURRENT MOD	SMA connector. Control input for analogue diode current modulation. $\pm 10\text{ V}$ input max, 40 mA/V sensitivity (modulation depth of 400 mA) and bandwidth of 40 kHz .
LASER DIODE	Female DE15 port to laser head. Provides diode current, I2C, 5 V and other control lines to the laser headboard.
INTERLOCK	The LDD will not power on the laser unless the pins on this port are shorted . A standard 3.5 mm audio plug is provided; see appendix D for pinout. The pins should be shorted using a relay for integration with laboratory interlock system; do not apply a voltage directly across the pins .
TEMPERATURE	Female DE9 port to laser head. TEC current and temperature sensor.
Ethernet	RJ45 10/100 twisted pair ethernet port.
USB	USB-A hi-speed (USB2) port.

1.3 Internal DIP switches

Special functions can be accessed by changing the state of dip switches within the unit. Note that these settings take effect on boot.

DIP 1	Access firmware update mode. Normally not required; Please contact MOGLabs if you have problems with updating firmware
DIP 2	Internal use only and should be left in the on position
DIP 3	Internal use only and should be left in the on position
DIP 4	The firmware will be restored to the original factory version using a stored image. Saved settings (stored in EEPROM) will also be restored to factory values: please check maximum allowed laser diode current, I_{lim} .

2. Menus

The MOGLabs LDD can be controlled via a detailed on-screen menu system. There are four push buttons to the right of the LCD display, to step through the menus. The large encoder knob can be rotated to change a value, and also acts as a push-button to change which character or digit is being changed.

The LCD display provides status and control information in the following pages:

Main	Critical values and settings
Settings	Full settings menu system
TEC	Temperature-related values and settings
Laser	Laser values and settings
Ramp	Piezo and ramp settings
Interlock	Interlock and photodiode settings
Ethernet	Ethernet (TCP/IP) configuration
Display	Display settings
System	System settings and debug information

2.1 Main menu

Laser Diode Driver	
Current Setpoint:	0.500 A
Current Limit:	0.500 A
Temp Setpoint:	+24.00 C
Temp Reading:	24.00 C
Settings ▶	
10.1.1.139:7802	24.00 C

The default screen on power-up, showing diode set current, current limit, set temperature, actual temperature. The color box shows that

the LD current is selected for adjustment. Rotate the adjust knob to change the current; use the ↓ button to select a different parameter, or to move the cursor down to the Settings→ line and press OK to enter into other menus.

2.2 Settings menu

Device Settings		Device Settings	
TEC ▶		Ramp ▶	
Laser ▶		Interlock ▶	
Ramp ▶		Ethernet ▶	
Interlock ▶		Display ▶	
Ethernet ▶		System ▶	
10.1.1.166:7802	20.00 C	10.1.1.166:7802	20.61 C

2.3 TEC menu

TEC	
Tec V:	0.00 V
Tec I:	0.00 A
PSU:	14.91 V
PSU Power:	0 W
Settings ▶	
10.1.1.146:7802	23.63 C

The TEC menu provides access to key operating parameters of the temperature controller, and a sub-menu to change TEC settings.

2.3.1 TEC settings

TEC Settings		TEC Settings	
Max I:	2.50 A	Temp Min:	+10 C
TEC Auto Polarity:	ON	Temp Max:	+35 C
TEC Invert:	ON	Kp:	125
Temp Min:	+10 C	Ki:	500
Temp Max:	+35 C	Kd:	400
10.1.1.166:7802	22.69 C	10.1.1.166:7802	22.69 C

PID coefficients can be adjusted here if needed. *Max I* sets an upper limit to the current to the TEC. *Temp Min* and *Temp Max* are limits; if the temperature falls outside this range, the temperature controller is disabled and current to the TEC is set to zero.

By default, the LDD will auto-detect the TEC polarity. Auto-polarity can be disabled and the polarity manually inverted with the *TEC Invert* menu option.

2.4 Laser menu

Laser Settings	
Setpoint I:	0.100 A
Measured I:	0.100 A
Measured V:	-0.17 V
PSU:	-5.63 V
PSU Power:	0 W
10.1.1.139:7802	24.03 C

Setpoint I displays the current requested when the LDD current output is turned on. *Measured I* and *Measured V* display the diode current and voltage measured at the LDD unit. *PSU* and *PSU Power* give indicative internal operating voltage and power of the current driver.

2.5 Ramp settings

Ramp Settings		Ramp Settings	
Ramp span:	0.00	Current ramp:	OFF
Ramp offset:	-1.00	Set current ramp ▶	
Ramp frequency:	9.0 Hz	Piezo enabled:	OFF
Ramp bias:	0 mA	Piezo mode:	RAMP
Current ramp:	OFF	Set piezo mode ▶	
10.1.1.139:7802	24.02 C	10.1.1.139:7802	24.04 C
Current Mod Select		Piezo Mode	
None		Internal Ramp	
External		External Piezo Mod	
Internal Ramp			
Negative Ramp			
10.1.1.139:7802	23.78 C	10.1.1.139:7802	23.59 C

The LDD has an internal ramp generator that can be used to drive a piezo and/or current ramp.

Ramp span, *Ramp offset* and *Ramp frequency* allow control of the ramp shape.

Current ramp indicates the *Current Mod Select* status selected from the *Set current ramp* submenu; *None*, *External*, *Internal Ramp*, *Negative Ramp*. If the current ramp is set to *Internal ramp* or *Negative ramp*, *Ramp bias* refers to the amplitude of the current ramp. An *External* current ramp can be provided via the rear SMA connector CURRENT MOD. Note that *Ramp bias* must be non-zero for the external CURRENT MOD to work.

Piezo enabled indicates *ON* or *OFF* status for the piezo signal passing through to the laser head. *Piezo mode* indicates *Set piezo mode* status; *RAMP* for *Internal Ramp* or *EXT* for *External Piezo Mod*. An external piezo modulation can be provided via the rear SMA connector PIEZO MOD.

2.6 Interlock

Interlock Settings		Interlock Mode	
Interlock State:	OK	Off	
PD Reading:	0.01 V	Case Switch	
PD Setpoint:	0.00 V	Photodiode	
Interlock Current:	1.000 A		
Interlock Mode:	OFF		
10.1.1.45:7802	20.00 C	10.1.1.139:7802	23.44 C

The LDD has a headboard-based interlock in the form of a case-switch and a monitor photodiode (see sec. 3.2). The interlocks are active only when the output current is greater than the *Interlock Current*.


When *Interlock Mode* is set to *Photodiode*, the interlock will cut the LDD current if *PD Reading* falls below *PD Setpoint*. To calibrate the *PD Reading* refer to section 6.2 of the MOGLabs MSA/MOA manual. *PD Reading* saturates at 4.94 V.

2.7 Ethernet settings

Ethernet Settings		Ethernet Settings	
Current IP:	10.1.1.146	IP Mask:	255.255.255.0
Static IP:	10.1.1.140	Gateway:	10.1.1.1
IP Mask:	255.255.255.0	Port:	7802
Gateway:	10.1.1.1	DHCP:	ON
10.1.1.146:7802	23.60 C	Restart ethernet	
		10.1.1.146:7802	23.65 C

TCP/IP network parameters. Note for these settings to take effect select the *Restart ethernet* menu option which is initially off-screen at the bottom.

2.8 Display



The LDD display settings, such as contrast and sleep timer can be adjusted here.

2.9 System information

System	System
Board Temps ▶	System Powers ▶
System Powers ▶	Board Voltages ▶
Board Voltages ▶	Versions ▶
Versions ▶	Restore Default Settings ▶
Restore Default Settings ▶	Firmware Update Mode
10.1.1.139:7802 24.07 C	10.1.1.139:7802 23.45 C
Board Temps	Sys Powers
PCB Temp: 36.75 C	TEC Power: 15 W
Tec CTRL: 34.19 C	Laser Power: 0 W
Laser CTRL: 34.56 C	Piezo Power: 0 W
Fan Speed: 0.0 %	Total Power: 15 W
10.1.1.139:7802 23.16 C	10.1.1.139:7802 24.17 C

Board Voltages		Version Information	
P3v3:	3.28 V	Serial number:	A07067
3v3 E:	0.00 V	Uptime:	631.64
P5v:	5.01 V	Microcontroller:	1.3.4
P14v:	14.05 V	Bootloader:	1.1.0
N14v:	-14.15 V	Display:	1.0.24
10.1.1.139:7802	23.79 C	10.1.1.139:7802	23.63 C

The System menu allows access to infrequently required information, including internal temperature readings and supply voltages from the main board as well as firmware versions. There are also options to restore default settings, Factory reset the settings and firmware, as well as restart the device.

3. Operation

3.1 Getting started

The MOGLabs LDD should be connected to the laser diode and TEC using heavy gauge screened cables, one for each of the current and temperature functions. The laser chassis connected to the controller should be in good electrical contact to an optical table or other surface, which in turn should be connected to the mains power supply electrical ground.

To operate the LDD:

1. Ensure the rear power switch is on, and the STANDBY/RUN switch is in the STANDBY position. In this mode, the current and temperature controller are disabled.
2. Switch from STANDBY to RUN. The TEC LED should turn yellow indicating that the temperature controller can now be enabled. If the indicator is still off, this indicates that an interlock or connection to the laser is missing.
3. If the key is switched back to STANDBY, both laser current and temperature controllers will be disabled.
4. Adjust the temperature setpoint as required, using the display or `mogldd` application. If the LDD was purchased as part of a laser system, this will already have been factory-set as required.

To use the display, move the cursor down to *Temp Setpoint* using the buttons on the right-hand side, and then turn the rotary encoder (ADJUST) to adjust the value.

5. Start the temperature controller by pressing the ON/OFF button below the TEMP LED. The LED below TEMP may briefly turn blue indicating that the the temperature controller is testing the polarity of the TEC. This process takes approximately three

seconds. Once the temperature controller is running the LED will turn green.

6. Adjust the current setpoint as required. It may be necessary to change the *Current Limit* to permit the desired setpoint current.
7. Switch the laser on by pressing the ON/OFF button below the CURRENT LED. The LED below CURRENT will turn blue, indicating that the current supply is starting. After three seconds the output current will ramp up to the desired setpoint value.
8. Once the current has hit the setpoint the LED will turn green, indicating the laser is operational.

Note that the temperature controller must be running, all interlocks must be enabled, and no errors can be present for the laser current to be enabled. Consult Appendix B for assistance with diagnosing issues.

3.2 Interlocks

The LDD includes a number of safety features for the protection of both personnel and equipment. These interlocks must be active during normal operation, as even a momentary interruption will trigger the safety procedure and disable the laser output.

The laser headboard interlocks (case switch and photodiode) only activate above an *interlock threshold current*, allowing for diagnostics and laser tuning to be carried out at lower currents (and hence lower power). These settings can be adjusted in the *Settings* → *Interlock* menu on the device, or through the headboard commands (section E.4).

3.2.1 Rear-panel interlock

The back-panel interlock connector is a 3.5mm audio plug whose pins must be shorted to enable the LDD. Intended for integration

with controlled-entry systems such as door interlocks.

3.2.2 Key switch

The key switch on the front-panel must be set to *Run* before the LDD output can be enabled. As a safety procedure, it is necessary to manually toggle the key switch after power-cycling the unit; it cannot simply be left in the ON position.

For integration into control systems, toggling the key switch can be performed through the command `TOGOVERRIDE`, however to be compliant with safety regulations this function should only ever be used in accordance with a user prompt confirming manual override.

3.2.3 Case interlock

High-power laser amplifiers produce a lot of scattered light, and are generally considered unsafe to operate with the lid open. The laser headboard includes a case switch that prevents the laser from being run above the interlock threshold current when the case is open.

3.2.4 Photodiode interlock

Some laser amplifiers include a photodiode sensor that monitors the approximate laser output power. An unexpected decrease in output power at a given current is an indication that the amplifier is operating unseeded, which can irreversibly damage the amplifier at high drive currents.

A threshold photodiode voltage can be specified (e.g. corresponding to 70% of the typical output power), below which the interlock will be tripped, helping to prevent damage from occurring. This setpoint can be specified using the `HB,VSET` command, and the currently-measured photodiode voltage can be queried using `HB,PDV`.

4. MOGLDD application

The host software program `mogldd` provides a graphical user interface that allows remote control of the LDD laser diode driver.

It may be necessary to install a firmware update (section 4.4) before being able to use the `mogldd` software. If the software detects an incompatibility it will offer to install the update, which can be obtained from the MOGLabs website.

4.1 Device discovery

Upon starting `mogldd`, a device discoverer is displayed (Figure 4.1) which searches for LDD devices. Simply select the device with the correct serial number and click *Connect*.





Figure 4.1: Example of the *device discoverer* window.

If the network does not permit device discovery and/or your LDD does not appear in the list, type the IP address of the unit (as displayed on the LCD screen) in the *Device address* box and click *Connect*.

4.2 Main display



The main window is divided into sections displaying the state of different aspects of the device, as follows. Values are periodically read from the unit and updated.

Power control Two buttons allow control of the TEC and laser diode current and show their status.

TEC Activates the TEC output and temperature feedback control.

Diode current Activates the laser diode output with soft start: the current increases slowly after a short delay.

Status Status of the front panel keyswitch and interlock, which must both be enabled to initiate the current controller.

Diode current

I_{actual} Measured current through laser diode.

I_{set} Desired output value of the current controller.

V_{diode} Measured voltage across laser diode.

Temperature The measured temperature of the laser, and the current and voltage of the TEC.

T_{actual} Measured temperature determined from the temperature sensor (either a thermistor or AD590 sensor).

T_{set} The setpoint temperature of the TEC controller.

I_{TEC} Measured current through TEC.

V_{TEC} Measured voltage across TEC.

4.3 Menu options

4.3.1 Actions


Command This opens a new dialogue for interactive command access to the device (Appendix E).

Restart LDD Resets the LDD device. Communications will be suspended. `mog1dd` will attempt to reconnect to the device for up to 10 seconds.

Restore settings Initiates restore of factory settings from EEPROM on the LDD device.

4.3.2 Config

Current Options for configuring the current controller, including features to prevent damage to attached devices.



The screenshot shows a dialog box titled 'Current' with a blue icon on the left and '?' and 'X' buttons on the right. It contains four rows of configuration options, each with a label, a text input field, and a unit. The first three rows have input fields with values: 'I set:' with '0.320', 'I limit:' with '2.20', and 'Maximum diode voltage:' with '5.0'. The fourth row, 'Cover interlock:', has a checked checkbox and the label 'Enabled'. The fifth row, 'Current modulation:', has an unchecked checkbox and the label 'Enabled'. An 'Exit' button is located at the bottom right of the dialog.

I set:	0.320	A
I limit:	2.20	A
Maximum diode voltage:	5.0	V
Cover interlock:	<input checked="" type="checkbox"/> Enabled	
Current modulation:	<input type="checkbox"/> Enabled	

I set An alternative control for setting the desired laser diode current.

I limit Defines an upper limit for the entered diode current on the LDD device.

- Max voltage** Defines upper limit on the diode voltage on the LDD device.
- Current mod** Enables external analogue modulation of the laser diode current via the rear panel modulation input.
- Temperature** The temperature settings dialogue allows setting critical temperature controller parameters.

T set An alternative control for setting the desired temperature.

Min/Max temperature, Max current Define the limits on the permitted setpoint temperature and the maximum TEC current. Exceeding these limits is indicative of equipment failure and will disable the device as a safety feature.

PID coefficients The LDD implements the feedback control via a standard PID (proportional integral differential) function:


$$u(t) = k_p e(t) + k_i \int_0^t e(\tau) d\tau + k_d \frac{de}{dt}.$$

$e(t)$ is the input error signal and $u(t)$ is the feedback response. k_p, k_i, k_d are scaled positive proportional (P), integral (I) and differential (D) gain coefficients. Typical values are $P = 125$, $I = 500$ and $D = 400$. New coefficients can be entered but are not activated until the adjacent Set button is pressed.

Inverted current Allows application of the LDD with a reverse polarity TEC.

Automatic inversion Instructs the LDD to determine the TEC polarity by checking which polarity causes a temperature decrease.

Network Allows configuration of network connection settings (IP address, mask, gateway and port). Particularly useful for configuring the network settings over USB. Note that changing the *Static IP* only has an effect if DHCP is disabled, or if DHCP name resolution fails. Note



The screenshot shows a window titled "Ethernet settings" with a close button (X) in the top right corner. The window contains a section titled "Configure ethernet settings". Below this title are several input fields and a checkbox:

- MAC address: 70:B3:D5:84:A2:06
- Current IP: 10.1.1.123
- DHCP: ☐ Enabled
- Static IP: 10.1.1.123
- Mask: 255.255.255.0
- Gateway: 10.1.1.1
- Port: 7802

At the bottom of the window are two buttons: "Save" and "Close".

that changing the ethernet settings will require the application to be restarted, and may also require the device to be rebooted. The port should be unchanged at 7802 to ensure that the `mogldd` suite of programs can continue to communicate with the device.

Reconnect Initiates device reconnection following communication interruption.

4.3.3 Help

About Displays device information.

Update Initiates firmware update (section 4.4).

Manual Opens web browser pointed at the MOGLabs support website to read the most up-to-date version of the manual.

4.4 Firmware update

From time to time, MOGLabs will release a firmware update that improves the device's functionality. The update is available from the MOGLabs website and should be installed on the device using the `mog1dd` application.

The firmware update procedure is started by selecting *Update* from the *Help* menu, or will be automatically activated if the software detects incompatibility with the current firmware version (Figure 4.2).




Figure 4.2: When `mog1dd` detects a version incompatibility it will offer to install a firmware update. Alternatively the *Update* option can be selected from the *Help* menu.

The update process depends on the hardware revision of the unit. Only Rev 4+ units are described by this manual, which are easily identified by the presence of a colour screen. A separate manual describes legacy units, which is available from the website.

It is strongly recommended to update the firmware using Ethernet, but it is possible to update over USB. When using USB it may be necessary to unplug the USB cable when the process appears to have become stalled at the “Waiting for reboot...” stage.

The LDD should not be in use while applying a firmware update; neither the front-panel interface nor a separate device connection should be used to interact with the unit simultaneously. The LDD must not be powered off or interrupted during the firmware update or the firmware could become corrupted.

Firmware for Rev 4+ units is distributed as a ZIP file that contains different firmware components. Upon opening the firmware update

tool (Figure 4.3), click the *Select* button and browse for this file. The tool will identify the components that need upgrading, which will be installed by clicking the *Update all* button.




Figure 4.3: The firmware update tool shows the current versions of the firmware components, indicating which are up-to-date (green) and which require updating (yellow or red).

The firmware update process typically requires the device to reboot several times. A prompt is displayed once the process is complete. Closing the firmware update tool will then allow `mog1dd` to be used.

A. Specifications

Parameter	Specification
-----------	---------------

Current regulator	
Output current	0 to 8 A
Max diode voltage	5 V up to 6 A
Display resolution	± 1 mA
Noise	TBD (10 Hz – 1 MHz)
Stability	Warmup time: 15 minutes
Analogue mod	40 mA/V, ± 10 V max, 40 kHz bandwidth

Temperature controller	
TEC current max	± 5 A
TEC voltage max	± 12 V
TEC power max	60 W
Stability	± 10 mK/ $^{\circ}\text{C}$
Sensor	NTC 10 k Ω
Range	-20 to $+80^{\circ}\text{C}$ standard
Display resolution	$\pm 0.01^{\circ}$

Note: The TEC is controlled with a linear regulator, which may overheat if the current load is high and the TEC voltage is low. Choose a TEC with 2 to 4 Ω resistance for optimal performance.

Parameter	Specification
-----------	---------------

Piezo driver	
Output	0 to 150 V, 10 mA
Bandwidth	1 kHz
Control	External analogue 12.5 V/V Internal 12-bit DAC (4096 steps)

Computer interface	
Ethernet	10/100 TP, RJ45
USB	USB2, plug type USB-A
Display type	2.8" 240x320 TFT

Connectors	
Current	DE15 high density 15-pin female
Temperature	DE9 9-pin female
Piezo	LEMO EPL.0S.302.HLN

Protection	
External interlock	3.5 mm audio plug (provided)
Cover interlock	Laser head cover interlock
Key interlock	STANDBY/RUN
Delayed soft-start	3 s delay + 1 s/A ramp (to 6A)
Open circuit	Laser cable, TEC, temperature sensor
Diode current	Digital setpoint limit I_{lim}
Photodetector	Seed failsafe

Parameter	Specification
-----------	---------------

Mechanical & power	
Fan	Dual 12 V DC ball-bearing
IEC input	80 to 264 47-63Hz
Dimensions	WxHxD = 250 × 79 × 292 mm
Weight	2.4 kg (excluding cables, laser head board)
Power	18 W (standby) 30 W to 100 W (low/high diode/TEC load)

B. Troubleshooting

The MOGLabs LDD detects a wide range of fault conditions and deactivates related circuitry accordingly. The front-panel LEDs provide indication of the state of these functions.

B.1 TEMP indicator

Colour	Status
DARK	Deactivated; check interlock and key switch
YELLOW	Standby mode, ready to be enabled
BLUE	Testing TEC direction
GREEN	Operational
RED	Error state; check screen for more information

B.2 CURRENT indicator

Colour	Status
DARK	Deactivated; check TEC controller
YELLOW	Standby mode, ready to be enabled
BLUE	Current ramping up to setpoint
GREEN	Current stable at setpoint
RED	Error state; check screen for more information

B.3 Error states

If an error is detected, the LCD display will provide an explanatory message from the table below.

Error Message	Solution
Internal PCB too hot	The internal circuit board has overheated. Check ventilation and ensure fan is not blocked.
Laser driver too hot	
TEC driver too hot	
Error reading temperature sensor	Unable to measure laser temperature. Check connections between LDD and laser headboard.
Laser temperature too high	TEC could not regulate temperature. Indicates incorrect TEC polarity, insufficient chassis cooling, or electrical failure of the TEC.
Laser temperature too low	
TEC error	Error controlling TEC
TEC open circuit	The TEC is not correctly connected, or has failed
TEC short circuit	
TEC max current too high	The voltage required to achieve the TEC current is too high, indicating the max TEC current must be reduced
Laser voltage too high	The voltage supplied to the laser diode has exceeded the maximum permitted voltage
Current control error	The current driver could not supply the required current.
Laser lid open	The laser lid has been opened, triggering the interlock (section 3.2)
Optical power too low	The photodiode is below threshold, triggering the interlock (section 3.2).

Error Message	Solution
Laser open circuit	The laser diode is open circuit. Check laser cable is correctly connected.
No interlock	The rear interlock is open circuit. Make sure the rear-panel interlock is shorted.
Remote interlock removed	Rear interlock was removed during operation.
Activate key first	Key needs to be in the RUN position.
Toggle key switch	The key needs to be turned off and on again before operation.
Key switch disabled	Key switch was turned off during operation.
Laser short circuit	The laser diode is short circuit. Check laser cable and diode.

For additional assistance please contact MOGLabs support describing the problem. Please ensure to include the device serial number and firmware versions.

C. Laser head board

A laser head interface board provides connection breakout to the laser diode, TEC, temperature sensor, and laser head interlock. It also includes a laser diode protection relay, a passive protection filter and a laser-on LED indicator.

A laser can be connected to the LDD controller via a B1048, B1055 or B1056 headboard (figs. C.1, C.2) which provides connections to the TA diode, TEC, and a passive NTC thermistor temperature sensor. All headboards have current capacity of 8 A for the laser diode and 5 A for the TEC. Connections from the B1048, B1055 boards to the components use Hirose DF59 “swing-lock” wire-to-board connectors; B1056 uses JST EH-series connectors. There is also a photodiode amplifier which can be used to generate a signal for monitoring the laser output. The boards include a solid-state protection relay, passive protection filters, and a laser-on LED indicator.

C.1 B1055/B1056 headboard

The B1055/B1056 are used in more recent MOGLabs amplifiers. They include connection to a cover interlock (pins must be shorted to enable laser diode current output). The switch SW selects cover interlock or photodiode interlock operation or both. Contact MOGLabs for further details if required, or refer to the MOGLabsMOA/MSA manual.




Figure C.1: MOGLabs MOA/MSA B1048 (superseded) laser head board showing Hirose DF59 “swing-lock” wire-to-board connector locations for diode (LD), NTC thermistor temperature sensor (Therm), TEC, photodiode, LED, and temperature/current head cables. The two TEC connections each have $+/-$ polarities and should be connected in parallel to provide greater current-carrying capacity. Each of the laser anode (+) and cathode (−) connectors have two pins connected in parallel. Each pin is rated to 2A.




Figure C.2: MOGLabs MOA/MSA B1055 (left) and B1056 (right) laser head board showing connector locations for tapered amplifier diode, temperature sensor, TEC, photodiode, cover interlock, LED and head cables.




Figure C.3: MOGLabs laser head board schematic (B1055).

D. Connector pinouts

WARNING: The CURRENT and TEMP connectors are intended for connection to a MOGLabs laser head board. They can provide high currents that may damage other devices. Considerable care should be taken if connecting to non-MOGLabs equipment.

Use only high quality cables (maximum length 3 m) with 26 AWG or larger conductors for current loads up to 4 A. Above 4 A, 22 AWG may be required (ask for advice from MOGLabs).

D.1 Current

Pin	Signal	Pin	Signal
1	Relay -	9	N/C
2	I2C SDA	10	Diode anode
3	Photodiode status	11	5v
4	Diode cathode	12	I2C SCL
5	Diode anode	13	Diode cathode
6	GND	14	Diode cathode
7	GND	15	Diode anode
8	GND		




Figure D.1: DE15, female, current supply connector (P6/P7) pinout. Relay (—) should be grounded to open the protection relay and enable operation. Pinout is the same on LDD and headboards. +5V is output on LDD, input on headboards.

D.2 Temperature

Pin	Signal
1	Thermistor +
2	GND
3	TEC -
4	TEC +
5	TEC +
6	Thermistor -
7	TEC -
8	TEC -
9	TEC +





Figure D.2: DE9, female, temperature supply connector (P1) pinout. Pinout is the same on LDD and on headboards, except that pin 2 is not connected on headboards.

D.3 Interlock



1	+3.3 V via 5 k Ω
2	No contact
3	Ground

Figure D.3: 3.5mm plug for Interlock. Either mono or stereo plug can be used. In both cases the tip and shell should be used, with the central connector unused for stereo plugs. Tip (pin 3) and shell (pin 1) should be short-circuited to enable current to the laser diode.

Note: Do not apply a voltage across the interlock pins, or the LDD may be irreversibly damaged!

D.4 Piezo



1	0 – 150 V
2	Ground

Figure D.4: LEMO EPL.0S.302.HLN connector for high voltage piezo actuator driver. Plug type FFA.0S.302.

E. Command language

The LDD is controlled via a text-based communications protocol which allows it to be easily integrated into existing control systems.

Commands are formed by a comma-separated list of parameters, as listed below. The parameters shown in square brackets are optional, and commands that are called without parameters are treated as queries.

Commands respond with a string that begins with either “OK” or “ERR” to indicate whether it was successful. Queries (such as measured temperature) return a value with units.

In some cases, the *actual* value may be different from the *requested* value due to parameter limits.

Note that the control commands for Rev 4+ hardware have a different structure to the legacy firmware. It is recommended to use these new commands, even where older commands are still recognised for backwards-compatibility purposes.

E.1 General functions

INFO INFO

Report device identification, including running firmware version and serial number. Please include this information in all correspondence with technical support.

VER VER

Report currently running firmware versions.

HELP HELP[,cmd]

Returns a list of available commands and short descriptive messages. Commands that have subcommands are indicated with >>, and the subcommands can be listed by specifying the `cmd` parameter.

STATUS **STATUS**

Responds with a description of the most recent error that occurred (as listed in Appendix B.3) or “No error”.

CLEAR **CLEARERROR**

Clears the most recent error as displayed on the screen and as returned by the **STATUS** query.

REPORT **REPORT**

Returns a dictionary (comma-separated list of key-value pairs) of present device settings. The order of entries in the dictionary is subject to change at any time.

INTERLOCK **INTERLOCK**

Report remote interlock status. Returns either ON indicating the interlock is present, or OFF when the interlock is missing.

KEY **KEY**

Report keyswitch status. Responds ON if the keyswitch is set to RUN, OFF if the keyswitch is set to STANDBY, or TOGGLE if the keyswitch must be manually toggled in accordance with safety guidelines.

TOGOVER **TOGOVERRIDE**

Overrides the requirement to toggle the keyswitch. In order to remain compliant with safety regulations, this command should only ever be issued in response to direct user intervention, such as a interactive alert or prompt.

BKLIGHT **BKLIGHT**[,value]

Screen back-light auto off time. **value** can be any number of seconds from 0 to 255. The back-light will remain on indefinitely with a setting of 0.

BTEMP **BTEMP**

Report measured board temperature.

FAN **FAN**

Report current fan speed, as a percentage of maximum.

E.2 Temperature control settings

TEC,ONOFF `TEC,ONOFF[,onoff]`

Activates the temperature controller (TEC) if `onoff` is ON, deactivates it if `onoff` is OFF, and reports the current status if `onoff` is not specified. Returns the string "AUTOPOLARITY" when performing the autopolarity test.

The keyswitch must be ON before the controller can be activated. Disabling the temperature controller will automatically disable the current supply.

TEC,TEMP `TEC,TEMP`

Reports laser temperature as measured by the headboard thermistor.

TEC,PID `TEC,PID,{P, I, D}[,value]`

Set or query the P, I or D gain constants used in regulating the temperature. These parameters should only need to be modified if the temperature sensor is physically relocated away from its factory position. `value` is rounded to the nearest whole number.

TEC,INVERT `TEC,INVERT[,value]`

Set or query whether to invert current to TEC. If the polarity of the TEC is incorrect, the device will heat when it is supposed to cool and the temperature will diverge. It is therefore important to ensure that current is inverted when the polarity is reversed from the default.

TEC,AUTOPOL `TEC,AUTOPOLARITY[,onoff]`

Enable/disable automatic determination of TEC polarity. This check takes approximately three seconds and is performed during the first activation of the temperature controller after power-up, as indicated by a blue LED on the front-panel.

TEC,TSET `TEC,TSET[,value]`

Set or query the desired temperature setpoint of the laser, as rounded to two decimal places. The optimal operating temperature is typically included as part of the laser diode specifications.

TEC,TMIN `TEC,TMIN[,value]`

Set the minimum permitted laser temperature, to prevent damage to the laser diode in case of control failure. `value` is a whole number. This temperature should never be reached if the TEC polarity is correct.

TEC,TMAX `TEC,TMAX[,value]`

Set the maximum permitted laser temperature, to prevent damage to the laser diode in case of control failure. `value` is a whole number. This temperature should never be reached if the TEC polarity is correct.

TEC,ILIM `TEC,ILIM[,value]`

Set the maximum TEC current or reports maximum TEC current if no value specified. A higher limit current may permit faster control convergence, but will require more heat to be dissipated through the baseplate.

TEC,ITEC `TEC,ITEC`

Reports the current measured through the TEC.

TEC,VTEC `TEC,VTEC`

Reports the voltage measured across the TEC.

TEC,THERM `TEC,THERM[,A,B,C][,val]`

Set or report the Steinhart-Hart coefficients for the temperature sense thermistor. For example, to set the coefficients for the standard 10 k Ω thermistor used by MOGLabs:

`TEC,THERM,A,1.124e-3`

`TEC,THERM,B,2.349e-4`

`TEC,THERM,C,8.49e-8`

To report one parameter, use for example `TEC,THERM,B`. To report all three use `TEC,THERM`.

E.3 Current control settings

CURRENT,ONOFF `CURRENT,ONOFF[,onoff]`

Activates the diode laser current supply if `onoff` is ON, deactivates it if `onoff` is OFF, and reports the status of the controller if `onoff` is not specified. The temperature controller must be running before the current supply can be activated.

CURRENT,ISET `CURRENT,ISET[,value]`

Set or query the setpoint current to be generated by the controller.

CURRENT,ILIM `CURRENT,ILIM[,value]`

Set or query the current limit. This limit is applied in hardware, and prevents damage to the laser diode by overdriving it.

CURRENT,MEAS `CURRENT,MEAS`

Reports the actual measured current through the laser. May differ from the setpoint current by a small amount.

CURRENT,MOD `CURRENT,MODULATION[,value]`

Set or query the current modulation mode. `value` can be one of the following options:

0, OFF, NONE No modulation

1, EXT, SMA External input (SMA connector)

+, INT Internal ramp (positive bias)

-, NEG Internal ramp (negative bias)

CURRENT,BIAS `CURRENT,BIAS[,value]`

Level of feedforward ("bias") applied to the current during the piezo sweep. `value` is a floating point value between 0 and 1.

E.4 Headboard settings

HB,HAS `HB,HAS`

Responds with whether the headboard is interlock-compatible.

HB,ISET *HB,ISET,[value]*

Threshold current above which the headboard interlock is activated. Below this current the headboard interlock is not active, allowing for diagnostics and adjustments to be carried out.

HB,INTER *HB,INTERLOCK[,mode]*

Set or query the headboard interlock mode. *mode* is one of the following strings

NONE Interlock disabled

CASE Case switch interlock

PD Photodiode interlock

If *mode* is not specified, the response includes the present interlock state, for example

- PD,1 indicates the photodiode interlock is active, and the interlock condition is currently being met (i.e. the measured optical power exceeds the threshold at this current).
- CASE,0 indicates the case interlock is in use, and the interlock condition is currently being violated (i.e. the case is open)

HB,PDV *HB,PDV*

Currently measured photodiode voltage (on compatible headboards).

HB,VSET *HB,VSET,[value]*

Minimum permitted measured photodiode voltage when drive current is above the threshold current.

E.5 Piezo driver settings

HV,ONOFF *HV,ONOFF[,onoff]*

Activates the high-voltage (piezo) driver. *onoff* can be ON or OFF; reports status of the piezo driver if parameter absent.

HV,SPAN `HV,SPAN[,value]`

Set or query the span setpoint. `value` will be rounded to one decimal place and should be between 0–120v.

HV,OFFSET `HV,OFFSET[,value]`

Set or query the high voltage offset. `value` can be any integer from –60 to 160.

HV,MOD `HV,MOD[,value]`

Set or query the piezo modulation mode. `value` is one of:

NONE No modulation

RAMP Internal ramp

EXT External input (via SMA connector)

HV,FREQ `HV,FREQ[,value]`

Piezo ramp frequency. `value` can be 1 to 100Hz; reports current value if parameter absent.

E.6 Ethernet settings

ETH,STATIC `ETH,STATIC[, "xxx.xxx.xxx.xxx"]`

Set IP default address based on decimal dotted-quad string (for example "10.1.1.180"). The double-quotes are required.

ETH,MASK `ETH,MASK[, "xxx.xxx.xxx.xxx"]`

Set IP mask based on dotted-quad string (for example "255.255.255.0").

ETH,GATE `ETH,GATE[, "xxx.xxx.xxx.xxx"]`

Set IP gateway based on dotted-quad string (for example "10.1.1.1").

ETH,MAC `ETH,MAC[, "xx:xx:xx:xx:xx:xx"]`

Set hardware MAC address based on colon-separated hexadecimal numbers (for example "70:b3:d5:84:a0:00").

ETH,PORT `ETH,PORT[,port]`

Set the TCP/IP port number for device communication.

ETH,DHCP **ETH,DHCP** [,onoff]

Enable or disable DHCP. Set to non-zero to enable DHCP; zero to use static IP address.

F. Communications

The LDD can be connected to a computer by USB or ethernet (TCP/IP) and integrated into existing control software. If you are experiencing difficulty in connecting to your LDD, please review the detailed instructions available at www.moglabs.com/ldd-software.html

F.1 Protocol

Communication follows a query/response protocol, where the user sends an ASCII string to the unit, and the unit sends an ASCII response to the originating source.

Messages are CRLF-terminated, so all statements must end with the carriage return (ASCII code 0x0D) followed by a new-line character (ASCII code 0x0A). Most terminal applications automatically append these characters. Responses from the unit should be buffered until CRLF is received. It is strongly recommended to check this response before continuing to ensure it does not indicate an error.

Statements are either **commands** or **queries**. A command is a statement that causes some action to occur, and the unit will respond with either "OK" or "ERR" depending on whether the command succeeded or not. For example,

```
> CURRENT,ON  
< ERR: TEC must be enabled first
```

```
> TEC,ON  
< OK
```

```
> CURRENT,ON  
< OK
```

Queries are statements to return a value, either returning the value requested or a message beginning with “ERR”. For example,

```
> TEMP
< 22.635 C

> TEMP
< ERR: Temperature sensor missing
```

F.2 TCP/IP

When ethernet is connected, the LDD will attempt to connect to the network using saved values, which can be altered using the front-panel menu system. If DHCP is enabled the device will first try to obtain an IP address via DHCP. If DHCP fails, the static settings will then be used. In both cases, the current address will be shown on the display.

F.2.1 Changing IP address

If your network does not use a DHCP server, you may need to manually change the IP address. This is easily done via the front-panel menu system in *Settings* → *Ethernet*.

F.3 USB

The LDD can be directly connected to a host computer using a USB cable (type A-male). The STM32 Virtual COM Port Driver device driver for the WindowsTM operating system is available from the MOGLabs website.

Connecting the LDD will install a new COM port on the machine. To determine the port number of the device, go to Device Manager (Start, then type `Device Manager` into the Search box). You should see a list of devices including “Ports” (Figure F.1).




Figure F.1: Screenshot of Device Manager, showing that the LDD can be communicated with using COM4. The port number might change when plugging into a different USB port, or after applying a firmware update.

The LDD can be identified as a COM port with the following name,
STMicroelectronics Virtual COM Port (COMxx)
where xx is a number (typically between 4 and 15).

If the port appears in Device Manager with a different name, then the driver was not successfully installed. If this occurs, disconnect the LDD from the host computer and reinstall the device driver.

